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### Instructions on environmental protection

Do not dispose of this equipment with normal household waste at the end of its life cycle; take it to a collection point for the recycling of electrical and electronic equipment.

It does not contain elements that are hazardous or toxic to humans but improper disposal will harm the environment.

The materials are recyclable as indicated on the marking. By recycling materials or other forms of reuse of old equipment, you are making an important contribution to the protection of the environment. Please contact your local authority for advice on collection points.

**8. Troubles and solutions**

No.	Problems	Possible Cause	Solutions
1	No ultrasonic activity	A. Power supply not connected B. Broken fuse C. Shorted cable D. Shorted transducer E. Broken integrated circuit board F. Other reasons	A. Check and plug power switch B. Check the fitted power supply and change the same specification fuse C. Connect fitted cable or replace a new D. Consult with your supplier E. Check the broken part and change it F. Consult with your supplier
2	Time control failure	A. Timer failure B. Other reasons	A. Replace the timer or digital control panel B. Consult with your supplies
3	No heating	A. Heater power switch had a bad linkage B. Fuse burn out C. Heating pad burn out D. Digital display board out of control E. Other reasons	A. Check heating plug and connect it well. Check outlet line with multimeter: OK if resistance value is a few hundred OHMs. B. Change the same specification fuse C. Replace the bad heating pad if open circuit D. Check the broken part and change it E. Consult with your supplier
4	Temperature control failure	A. Thermostat loosen B. Thermostat tube broken C. Digital display out of control D. Other reasons	A. Fasten the thermostat header B. Replace the thermostat C. Check the broken part and change it D. Consult with your supplier
5	Not well cleaning	A. Ultrasonic cleaning is not efficient B. Check liquid level in tank C. Temperature too high or too low D. Unsuitable cleaning fluid E. Other reasons	A. Adjusting the ultrasonic power with the control knob B. The recommended level is 2/3 of the tank capacity. Never less than 7 cm. C. Adjust the temperature to the most suitable one D. Stop and disconnect the power supply, replace the appropriate liquid after previous cooling of the liquid. E. Consult your supplier
6	Electric Leakage	Machine not grounded	Check if machine earth wire loosens

**INDEX OF LANGUAGES**

Spanish ..... 2-10  
 English ..... 11-19  
 French ..... 20-28

**TABLE OF CONTENTS**

Features ..... 11  
 1. Ultrasonic principle ..... 12  
 2. Preparation ..... 12  
 3. Safety precautions ..... 13  
 4. Applications ..... 14  
 5. Digital ultrasonic cleaner operation, structure and accessories ..... 14  
 6. Different ways of cleaning ..... 16  
 7. General specifications ..... 17  
 8. Troubles and solutions ..... 18

Thank you for purchasing the ultrasonic cleaner. Please take the time to read these operating instructions before use and retain them for future reference. Failure to follow these instructions may lead to serious personal injury and damage to property.

**Features**


- Stainless steel housing, tank and lid.
- Stamp tank without welding gap for better waterproof.
- Cooling fan
- Moisture-proofed PCB
- Industrial grade IC & ultrasonic transducer

**PLEASE READ THIS INSTRUCTION MANUAL CAREFULLY BEFORE USING THIS APPLIANCE TO GET THE BEST PERFORMANCE AND LONGEST LIFE FROM IT.**

## 1. Ultrasonic principle

Ultrasonic cleaning is based on the cavitation effect caused by high frequency ultrasonic wave vibration signal in the fluid. Microscopic bubbles are formed, and they implode violently causing the cavitation which create an intense scrubbing action on the surface of the item being cleaned. The bubbles are small enough to penetrate microscopic crevices, cleaning them thoroughly and consistently. Ultrasonic cleaning is extremely effective at removing dirt and grime which would normally require tedious manual cleaning by hand. It has been used to clean a wide variety of instruments and mechanical parts such as carburetors, returning them to almost "like new" condition without damage to delicate parts.

## 2. Preparation

- Carefully unpack the cleaner and remove all traces of packing materials from it. Visually inspect the cleaner for any parts that may have become loose or damaged during transit.
- Place the ultrasound on a clean, flat surface where the equipment receives adequate ventilation and ensure all controls are turned off. Also ensure the drain tap is in the closed position (in the references that include it).
- Ensure the power lead is securely plugged into the cleaner and that no part of the lead is likely to come into contact with moisture.
- Carefully fill at least 1/2 of the tank with a solvent solution. For demanding cleaning we recommend a small amount of washing up liquid, this will help increase the cleaning performance. The cleaner is now ready for use.

### Content

- |                 |                    |
|-----------------|--------------------|
| a: Main machine | b: Soundproof lid. |
| c: Power lead   | d: Mesh basket     |
| e: Manual       |                    |

While the machine is working normally, the ultrasound and tank tuning gives a well graded sound, and no tremors on the water surface. If there are discontinuous waves, please add or release some washing solution in the tank; stopping waves is better for cleaning objects.

On the condition of making sure the object is cleaned, please run the machine discontinuously as far as possible (no more than 30 minutes), long-time operation will cause the temperature of the casing.

## 6.2. For better cleaning effect

- The item should be immersed well into the water (do not exceed the "MAX" mark)
- Warm water (never hot or boiling) and a small amount of washing up liquid will help increase the cleaning performance.
- Make sure there is good space around each part you put in the tank, The more items that you place in the tank the less efficiently it will clean. It is not advised to overlap items, If lay items one on top of the other the system doesn't work the way it's supposed to.
- Basket use. Do not put items directly into the bottom of the tank, better cleaning effect with a suitable basket. A metal basket absorbs less ultrasonic energy at about 8%.
- Temperature: The higher temperature the better cleaning effect. However, when temperature exceed 70°C ~ 80°C, cleaning effect will be affect. The best temperature we suggest is 40°C ~ 60°C.

## 7. General specifications

Power supply:	220VAC 50Hz
Temperature range:	0-80°C (from ambient temperature to 80°C) 1°C in 1°C selection
Time range:	1 – 99 minutes; selection from 1 in 1 minute
Display	LED, two four-digit LEDs (red colour)

Code	Frequency (Hz)	Capacity (L)	Drainage tap	Tank size LxWxH (mm)	Unit size LxWxH (mm)	Heating power (W)	Ultrasonic power (W)	Timer (min)	Temp (°C)
ZCC016	40	2	No	150x135x100	175x160x210	100	60	1-99	RT-80
ZCC017	40	3.2	No	240x135x100	265x165x220	100	120	1-99	RT-80
ZCC019	40	6.5	Yes	300x150x150	325x180x280	300	180	1-99	RT-80
ZCC020	40	10	Yes	300x240x150	325x265x280	300	240	1-99	RT-80
ZCC021	40	15	Yes	330x300x150	360x325x285	400	360	1-99	RT-80
ZCC018	40	22	Yes	500X300X150	550X330X280	500	480	1-99	RT-80
ZCC022	40	30	Yes	500x300x200	530x325x325	500	600	1-99	RT-80

## 6. Different ways of cleaning

- **General Cleaning** — use only tap water. Use warm temperature about 50°C.
- **Enhanced Cleaning** — add few drops of standard washing up liquid, liquid soap, or detergent into water, other non-acidic cleaning agents can also be used.
- **Extensive Cleaning** — removing tarnish, carbon & rust from non-plated metals, it is recommended to use specialist cleaning solution in associated using ultrasonic cleaner.

**WARNING:** Strong acid or alkaline cleaning solution will cause corrosion, rust and even puncture of tank or machine body. To overcome this problem, dilute to mild PH solution or request for tank made of specific grade of stainless steel like SUS304.

The cleaning solution will deteriorate in effectiveness over time and use. It is important to regularly change the fluid and carefully wash the inside of the cleaning tank in order to preserve the effectiveness and longevity of the cleaner. Do not use corrosive or abrasive cleaning products on the tank, which must be wiped down and dried before it can be re-connected to the electrical supply.

**NOTE:** If the machine starts to spark, smoke, smell of burnt electrics or displays any other fault the operator must immediately stop the machine, isolate the electrical supply and contact the supplier. Continued use is dangerous.

### 6.1. Advantage

Ultrasonic is widely used throughout industries to remove difficult contaminants from the parts during or after manufacturing process which might require a stage of cleaing before the next process. In general, if an item can be cleaned with liquid, it can be cleaned much faster and more thoroughly with an ultrasonic cleaner, Compared to traditional solvent/scrubbing methods of cleaning in a parts washer, our ultrasonic cleaners:

- Are more effective at removing contamination
- Are quicker to achieve the results you want
- Save employee time (and consequently money) by enabling 'switch on and leave' cleaning
- Are heated to allow cleaning solvents to work more effectively
- Use mechanical timers for simplicity, precision and longevity
- Produce less mess during use
- Less chemistry
- Reproducible result

## 3. Safety precautions

**Keep it away from children!**

This device is not intended to use by individuals with restricted physical sensory or metal capacities or those with lack of experience or knowledge, include children, unless they are supervised by an individual who is responsible for their safety or have received training in operating the device.

**Please read the following very carefully as failure to comply may invalidate your guarantee**


1. DO NOT run the cleaner continuously for more than one hour at a time, as doing so can damage the internal components.
2. DO NOT operate the unit without fluid in the tank. Always ensure the fluid is no higher than the max mark. Always ensure there is a minimum depth of 7cm. Ideally 2/3 of the tank capacity should be filled.
3. Strong acids or flammable cleaners are strictly prohibited. If solvent type cleaner is added in the tank, please stop working for a while after the 30-minute cycle.
4. DO NOT drop any item into the tank as this may cause damage to the transducer. Always place the items gently into the tank and use the basket whenever possible.
5. The more items that you place in your cleaning bath the less efficient it will clean.
6. It is not advised to overlap items. Always allow plenty of clear space between the items.
7. Do keep cap on during use. This will prevent splashing and reduce liquid evaporation as well as reduce equipment noise.
8. Never immerse the machine or power cord in water or other liquid.
9. DO NOT touch the power plug with wet hands, especially when inserting or removing the plug.
10. DO NOT touch the unit if the machine has fallen into water during operation. Remove the power plug from the socket first.
11. DO NOT disassemble the machine, except by professionals.
12. UNPLUG the power source while filling or emptying the tank.
13. DO NOT spray water or liquid over the device and the control panel
14. DO NOT operate the cleaner without proper grounding.
15. DO NOT place the device on a soft surface, where the vents could be blocked.
16. Always turn the heater off after using as leaving it on can cause the fluid to evaporate and damage the internal components.
17. Upon completion of the cleaning cycle, turn the heater button off and isolate the machine from the electrical supply.
18. Take care when adding or removing items from the cleaning tank as the fluid is likely to be hot and displaced fluid can damage the internal components. Any displaced fluid must be dried up immediately.
19. In the event of failure/emergency, disconnect the mains supply by removing the plug from the mains socket.

## 4. Applications

This list is almost endless. Provided the product is non porous and can normally be immersed in water almost anything can be thoroughly cleaned. Here are some examples:

- Jewellery especially gold, silver & platinum
- Watchstraps
- Coins and other collectibles
- PCB Boards etc
- Engine/Model parts
- Toothbrushes & Dentures
- Electrical components
- Diesel injection pumps
- Printer heads and toner cartridges
- Motorcycle radiators
- Vehicle differentials
- Milking parlour equipment
- Horse bits&stirrups &horse brasses
- Tattoo needles
- Sugical equipment
- Motorcycle engine crank cases
- Engine cylinder heads
- Turbochargers
- Bicycle derailleurs
- Ultrasonic cleaning is not recommended for the following gemstones: Opal, Pearl, Emerald, Tanzanite, Malachite, Turquoise, Lapis and Coral.

## 5. Digital ultrasonic cleaner operation, cleaner structure and accessories

### 1. Fill stainless steel tank with liquid

Plug the ultrasonics into grounded outlet; After powing on, The temperature displays the actual environment temperature, LED displays 3 minutes, which is default time.

### 2. Turn on the power switch

Plug the ultrasonics into a grounded outlet. Turn on the power switch on the back of the tank. ("O" is closed, and "I" is open)

### 3. Start ultrasonic cleaning

For the first boot, time display show setting time at random, press time key+ one time, "the setting time" on the display will be increase by 1 minute correspondingly, touch time key- to decrease time setting by 1 minute. Touch ON/OFF switch after you have properly set the time, ultrasonic indicator will turn red and the machine will start working based on the parameters set. Meanwhile the timer will begin to count down and stop working when the time shows "00:00", automatically.

### 4. Time setting:

Press the "TIME-/TIME+" shortly at a time means time increase/reduce 1min at a time, press and hold means time will increase/reduce by 1mins continuously. The cleaning stops when the timer counts down to 00:00. If the unit needs to be stopped before the timer runs up, press ON/OFF button.

**NOTE: Do not run the cleaner continuously for more than one hour**

### 5. Set temperature:

Press "temperature-/temperature+" shortly at a time means time increase/reduce 1°C at a time, press and hold there means time will increase/reduce by 1°C. Usually, the best results are within 40°C ~ to 60°C.

6. After time and heat setting, press ON/OFF key to start both function. Press "Heating" button "OFF", the heating function will stop. The same for Ultrasonic function.

7. Empty the tank and clean both the outside and inside of the cleaner with clean and dry cloth for next use.

**NOTE: Do not pour water out until it's cooling, Hot water will hurt you and damage the machine itself.**