

Elevadores de laboratorio | Jacks

Marca|Trademark:**Nahita**

Diseñados para elevar y mantener todo tipo de equipos y aparatos como frascos, cubetas, baños, agitadores, placas calefactoras, etc. Fabricados en aluminio y con soporte tipo tijera, pueden ser regulados en altura mediante un mando manual situado en el frente del elevador y disponen de 4 patas de goma antideslizantes para asegurar una buena estabilidad y evitar movimientos bruscos. Los modelos de 15 y 20 cm de lado presentan un orificio en uno de los lados de la plataforma para introducir una varilla que sirva de soporte a aros, pinzas etc.

Nahita laboratory jacks are designed to elevate and support different equipments and items such as flasks, bottles, jars, baths, stirrers, heating plates, etc. Made of aluminium and with a scissor-type support they can be height adjusted by a knob located at the front of the apparatus and they are provided with 4 anti-skid legs to assure a good stability and avoid sudden movements. Models 63217015 and 63217020 present a small hole at one corner of the platform to accommodate a frame rod (not included) for mounting ring stands, clamps and others.



Referencia Code	Plataforma Platform	Altura Height	Carga máxima Max. load
63217010	100x100 mm	40-140 mm	14 Kg
63217015	150x150 mm	50-270 mm	27 Kg
63217020	200x200 mm	50-320 mm	38 Kg