



-45°C Chest Freezers

Model	LT C150	LT C300	LT C400	LT C500
Technical				
Gross volume:	140 litres	296 litres	383 litres	495 litres
Temp. range:	-10°C to -45°C	-10°C to -45°C	-10°C to -45°C	-10°C to -45°C
Refrigerant:	R404a	R404a	R404a	R404a
Baskets:	1	1	1	3
Power consumption 24/h:	3,5 kWh	5,3 kWh	5,8 kWh	6,5 kWh
Voltage:	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz
Compressor type:	Single stage	Single stage	Single stage	Single stage
Features				
- Audible and visual alarm	•	•	•	•
- High/low alarm	•	•	•	•
- Battery backup (optional)	•	•	•	•
- Prepared for CO2 backup	•	•	•	•
- Castors	•	•	•	Legs
- Handle lock	•	•	•	•
Outer dimensions				
Height:	885 mm	885 mm	885 mm	890 mm
Width:	720 mm	1260 mm	1560 mm	1660 mm
Depth:	695 mm	695 mm	695 mm	758 mm
Gross weight:	59 kg	82 kg	95 kg	122 kg
Net weight:	55 kg	76 kg	88 kg	107 kg
Product no.	N111002	N111003	N111004	N111005

Storage capacity	LT C150	LT C300	LT C400	LT C500
Box height	Total no. of boxes	Total no. of boxes	Total no. of boxes	Total no. of boxes
50 mm	84	216	282	315
75 mm	60	156	204	228
100 mm	39	99	129	144

For further information and details please refer to “Accessories”.