



-60°C Chest Freezers

Model	XLT C75	XLT C150	XLT C300	XLT C400	XLT C500
Technical					
Gross volume:	74 litres	140 litres	296 litres	383 litres	495 litres
Temp. range:	-30°C to -60°C	-30°C to -60°C	-30°C to -60°C	-30°C to -60°C	-30°C to -60°C
Refrigerant:	Nature R	Nature R	Nature R	Nature R	Nature R
Baskets:	1	1	3	4	3
Power consumption 24/h:	4,0 kWh	3,5 kWh	6,2 kWh	7,7 kWh	6,0 kWh
Voltage:	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz
Compressor type:	Single stage	Single stage	Single stage	Single stage	Single stage
Features					
- Audible and visual alarm	•	•	•	•	•
- High/low alarm	•	•	•	•	•
- Battery backup (optional)	•	•	•	•	•
- Prepared for CO2 backup	•	•	•	•	•
- Castors	•	•	•	•	Legs
- Handle lock	•	•	•	•	•
Outer dimensions					
Height:	895 mm	885 mm	885 mm	885 mm	890 mm
Width:	550 mm	720 mm	1260 mm	1560 mm	1660 mm
Depth:	645 mm	695 mm	695 mm	695 mm	758 mm
Gross weight:	60 kg	59 kg	82 kg	95 kg	122 kg
Net weight:	55 kg	55 kg	76 kg	88 kg	107 kg
Product no.	N112001	N112002	N112003	N112004	N112005

Storage capacity	XLT C75	XLT C150	XLT C300	XLT C400	XLT C500
Box height	Total no. of boxes	Total no. of boxes	Total no. of boxes	Total no. of boxes	Total no. of boxes
50 mm	40	84	216	282	315
75 mm	25	60	156	204	228
100 mm	20	39	99	129	144

For further information and details please refer to "Accessories".