



-86°C Chest Freezers

Model	ULT C75	ULT C200	ULT C300	ULT C400
Technical				
Gross volume:	74 litres	198 litres	296 litres	383 litres
Temp. range:	-60°C to -86°C	-60°C to -86°C	-60°C to -86°C	-60°C to -86°C
Refrigerant:	Nature R	Nature R	Nature R	Nature R
Baskets:	1	2	3	4
Power consumption 24/h:	5,1 kWh	6,0 kWh	7,4 kWh	8,1 kWh
Voltage:	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz	220V-240V 50/60Hz
Compressor type:	Single stage	Single stage	Single stage	Single stage
Features				
- Audible and visual alarm	•	•	•	•
- High/low alarm	•	•	•	•
- Battery backup	•	•	•	•
- Prepared for CO2 backup	•	•	•	•
- Castors	•	•	•	•
- Handle lock	•	•	•	•
Outer dimensions				
Height:	895 mm	885 mm	885 mm	885 mm
Width:	550 mm	920 mm	1260 mm	1560 mm
Depth:	645 mm	695 mm	695 mm	695 mm
Gross weight:	60 kg	69 kg	82 kg	95 kg
Net weight:	55 kg	64 kg	76 kg	88 kg
Product no.	N113001	N113002	N113003	N113004

Storage capacity	ULT C75	ULT C200	ULT C300	ULT C400
Box height	Total no. of boxes	Total no. of boxes	Total no. of boxes	Total no. of boxes
50 mm	40	117	216	282
75 mm	25	84	156	204
100 mm	20	54	99	129

For further information and details please refer to "Accessories".