

# Upright freezer ULT U250

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable

## Specifications

### Technical

Gross volume:	253 litres
Temp. range:	-60°C to -86°C
Refrigerant:	Nature R
Inner doors:	2
Adjustable shelves:	5
Power consumption 24/h:	9,0 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage

### Features

- Temperature logger with USB port
- Door alarm
- High/low alarm
- 48h battery backup
- Castors

### Outer dimensions

Height:	2055 mm
Width:	600 mm
Depth:	595 mm
Gross weight:	135 kg
Net weight:	121 kg

### Inner dimensions

Height:	1420 mm
Width:	420 mm
Depth:	370 mm

Product no.	N123002
-------------	---------



## Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	4 x 2 = 8	3	24	N321200
	5 x 2 = 10	12	120	N321210
			<b>= 144</b>	
75 mm	2 x 2 = 4	3	12	N321300
	3 x 2 = 6	12	72	N321310
			<b>= 84</b>	
100 mm	2 x 2 = 4	3	12	N321400
	2 x 2 = 4	12	48	N321400
			<b>= 60</b>	