

# Rheomat R120

- Rotational viscometer for measuring viscosity according ISO and DIN.
- Rheomat R 120 excels in easiest handling and a large measuring range.
- Three fixed rotational speeds according application.
- Only one measuring system, therefore no errors due to wrong measuring system.
- Display of viscosity in Pas or mPas.

## Fields of operation

- Fast viscosity measurement
- Control of incoming goods
- Quality control
- And a lot more.

## Technical Data

Measuring principle	Rotational viscometer
Rotational speed range	20 to 1000 min <sup>-1</sup>
Torque range	up to 10,0 mNm
Viscosity range	20 mPas to 10 <sup>6</sup> Pas (Depending on measuring system)

Rheomat R 120 is delivered with its measuring system and power supply in a carrying case.

The Set includes:

- Viscometer R 120
- Power supply
- One proRheo measuring system
- Carrying case

