

Adapter



Adapter Thermometer Pocket

- Glass pocket for standard laboratory thermometers to be incorporated into glassware assemblies
- Closed end, made to fit stillheads

Item Code	Cone Size	Approx Stem Length (mm)	Pack of
042.202.01	14/23	45	10
042.202.02	19/26	45	10
042.202.03	24/29	45	10

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



Adapter, Air Leak Tube or Gas Inlet Tube

- Accurate dimensions
- compact and lightweight
- Durable
- Resistant to thermal shock
- Used to monitor and detect gas leaks in the apparatus

Item Code	Cone Size	Pack of
044.202.01	14/23	10
044.202.02	19/26	10
044.202.03	24/29	10
044.202.04	29/32	10

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



Adapter, Cone Screw Thread

- Ideal to securely connect thermometers and inlet tubes to glassware assemblies
- With a open top compression cap and O-ring

Item Code	Cone Size	Acceptable Dia through Screw Thread mm	Pack of
056.202.01	14/23	6	10
056.202.02	19/26	6	10
056.202.03	24/29	6	10
056.202.04	29/32	6	10

ASTM Sizes

Item Code	Cone Size	Pack of
056.470.01A	14/20	1
056.470.02A	19/22	1
056.470.03A	24/40	1

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



ENQUIRE