




NEW ROTOR  
4 X 100 ML  
Cat. N° 12200



**MPW**  
MED. INSTRUMENTS  
REF 13200  
SN 003/19  
001  
MAX 290g

	Type	Amount of tubes per NEW rotor 4x100	Amount of tubes per OLD rotor 4x100
	<b>Rotor 12200 4 x 100 ml</b>		
13113 	<b>Hanger 12 x Vacuette/Vacuteiner® 13 x 75 mm</b>	<b>48</b>	<b>48 Vacuette 32 Vacutainer</b>
13200 	<b>Bucket 100 ml</b>		
14013 	Carrier 8 x Vacuette/Vacuteiner® 13 x 75/100 mm	32	16
14016 	Carrier 7 x Vacuette/Vacuteiner® 16x100 mm	28	16
14020 	Carrier 5 x 10 ml Nalgene® tube	20	16
14021 	Carrier 10 x 2/1,5 ml Eppendorf® tubes	40	20
14023 	Carrier 1 x 30 ml Nalgene® tube	4	4
14026 	Carrier 1 x 50 ml skirted tube	4	4
14026C 	Carrier 1 x 50 ml Falcon® tube	4	4
14028 	Carrier 50 ml round bottom glass tube	4	4
14029 	Carrier 12 x 7ml tubes	48	20
14100 	Carrier 1 x 100 ml tube	4	4
13201C 	<b>Bucket 100 ml with cap</b>		
14013 	Carrier 8 x Vacuette/Vacuteiner® 13 x 75/100 mm	32	16
14016 	Carrier 7 x Vacuette/Vacuteiner® 16 x 100 mm	28	16
14020 	Carrier 5 x 10 ml Nalgene® tube	20	16
14021 	Carrier 10 x 2/1,5 ml Eppendorf® tubes	40	20
14023 	Carrier 1 x 30 ml Nalgene® tube	4	4
14026C 	Carrier 1 x 50 ml Nalgene® tube	4	4
14028 	Carrier 1 x 50 ml round bottom glass tube	4	4
14029 	Carrier 12 x 7ml tubes	48	20
14100 	Carrier 1 x 100 ml tube	4	4
13215 	<b>Hanger 2 x 15 ml Falcon® tubes</b>	<b>8</b>	<b>8</b>



**MPW MED. INSTRUMENTS**  
[www.mpw.pl](http://www.mpw.pl)

46 Boremlowska Street  
04-347 Warsaw Poland

Phone: +48 22 610 56 67  
Phone: +48 22 740 20 25  
Fax: +48 22 610 81 01  
E-mail: [mpw@mpw.pl](mailto:mpw@mpw.pl)



Technical specification subject to change. MPW MED. INSTRUMENTS is not responsible for printing mistakes. For further accessories please contact us directly or see at [www.mpw.pl](http://www.mpw.pl)