

SMZ140 Series

SMZ-143-N2LED





SMZ140 Series

SMZ-143-N2LED

Optical system	Greenough
Observation tube	Trinocular head
Inclination	45° inclined, 360° rotating
Trinocular light split	100:0/0:100 on right tube
Interpupillary distance	54-76mm
Diopter adjustment	On both tubes, +/- 5 diopter
Eyepieces	Widefield WF10X/20mm
Objectives system	Zoom, ratio 4:1
Objectives	1X~4X
Working distance	80mm
Stand type	Pole type stand with transmitted illumination
Head holder	For Ø25mm column and Ø74mm head with incident illumination
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stroke	45mm
Incident illumination	LED 3W with intensity control
Transmitted illumination	LED 3W with intensity control
Transformer	External
Power supply	100-240V (CE)
Accessories included	Dust cover, black/white and frosted stage plates
Dimensions LxWxH	237x170x335mm
Base	237x170mm
Column/Arm	Height 210mm, Ø25mm
Net weight	6.2kg