



Адва

AD110

**Standard Professional
pH-TEMP
Portable Meter**

AD110 is compact pH meter featuring automatic pH calibration, HOLD function, low battery indication, temperature measurement and compensation.

		AD110
Range	pH Temp	-2.00 to 16.00 pH -5.0 to 105.0 °C
Resolution	pH Temp	0.01 pH 0.1 °C
Accuracy	pH Temp	± 0.02 pH ± 0.5 °C up to 60 °C/ ± 1 °C outside
pH Calibration		automatic, 1 or 2 points, with 3 memorized buffers value
Temp. calibration		factory
Temp. compensation		automatic or 25 °C fixed
Probe		pH electrode with BNC (AD1230B) temp. probe (AD7662)
Battery type / life		4 x 1,5V AA / approx. 1500 hours of use
Auto-off		after 8 minutes of non-use
Environment		0 to 50 °C; RM max 95%
Dimensions / weight		70 x 145 x 26 mm / approx. 220