

Adwa



Professional Multi-Parameter Splash-Proof pH-ORP-TEMP Portable Meter

AD130 is a heavy-duty portable pH meter designed to provide laboratory results and accuracy even under harsh industrial conditions.

		AD130
Range	pH ORP Temp	-2.00 to 16.00 -1999 to +1999 mV -20.0 to 120.0°C (-4.0 to 248.0°F)
Resolution	pH ORP Temp	00.01 0.1 up to ±699.9 mV / 1 outside 0.1 °C (0.1°F)
Accuracy	pH ORP Temp	± 0.01 pH ±0.2 up to ±699.9 mV / ±1 outside ±0.4°C (±0.8°F)
Calibration		pH up to 2-points
Buffers		pH 1.68, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45
Temperature compensation		ATC/MTC (-20 to 120°C)
Calibration ORP		factory
UNIT		°C and °F
LCD		2 rows
Hold function		yes
Low battery indication		yes
Stability indicator		yes
Auto shutt-off		20 min. (can be disabled)
Power supply		4 x 1.5V AA batteries
Battery life		500 hours approx