



FLASK WITH DIVIDED  
NECK

KOLBEN MIT TEILUNG  
IM HALS

Cat. No. 1520

Catalog No. Katalog-Nr.	Capacity (ml) Kapazität	h (mm)	Packing (pcs) Packung (Stück)
1632431560725	50 + 5	155	10
1632431560830	100 + 10	205	10
1632431560837	200 + 20	245	10
1632431560938	250 + 25	280	10
1632431561143	500 + 50	320	8
1632431561144	1000 + 100	400	6
1632431561145	1000 + 10	300	7

ISO  
3507



BOTTLE, SPECIFIC  
GRAVITY, GAY-LUSSAC

PYKNOMETER NACH  
GAY-LUSSAC

Cat. No. 1621

Catalog No. Katalog-Nr.	Capacity (ml) Kapazität	Accuracy limits $\pm$ ml Genauigkeitsgrenze	d (mm)	h (mm)	SJ NS	Packing (pcs) Packung (Stück)
1632431722019	10	1.000	27	55	7/11	10
1632431722023	25	2.000	40	75	10/13	10
1632431722025	50	3.000	50	85	10/13	10
1632431722030	100	5.000	57	100	10/13	10

ISO  
3507

