

TECHNICAL SPECIFICATIONS

Trade name	ELECTRONIC TIMER
Use	Timing up to 19 hours 59 minutes 4 digit
Material	ABS case
Others	Electronic with alkaline batteries

DIMENSIONAL SPECIFICATIONS

Dimensions,	mm. 50,00 x 62,00 x 23,00
-------------	---------------------------



GENERAL FEATURES

CE IVD	No
--------	----

OTHER FEATURES

Instructions	Included
Other	Magnet for metal surface attachment
	Clips converts to stand

PACKAGING

Case	1 pc.
------	-------