

Myr

MCUT MICROTOME DISPOSABLE BLADES



MCUT microtome disposable blades are aimed to section the specimens at a certain selected stable thickness through required technical expertise such as sophisticated heating processes and polishing/coating techniques. The end result is exceptional sharpness, durability and consistency. MCUT, stainless steel microtome blades, allow users to deliver precise sections throughout all sample types.

Latest plasma technology with unique hardened cutting edge, improving the sharpness and durability of the blades and creating a superior cutting performance.

Available in two models, MCUT-3508 and MCUT-3514, low-and high profile universal blades that can be used with any microtome and cryostat in the market.

REF.	Description	Edge angle	length x height x thickness
MCUT-3508	Disposable microtome blades (low profile)	35°	80 x 8 x 0,254 mm
MCUT-3514	Disposable microtome blades (high profile)	35°	80 x 14 x 0,317 mm