

## WOULF BOTTLES WITH NECK AND OUTLET

Necks for connection to plant.


## WOULFFSCHE FLASCHEN mit Hals und Bodentubus

Cat. No. 2035						
Catalog No. Katalog-Nr.	Capacity (ml) Kapazität	d (mm)	h (mm)	SJ1	SJ2	Packing (pcs) Packung (Stück)
1632415041500	500	95	145	19/26	19/26	4
1632415041940	1000	110	180	24/29	19/26	4
1632415041950	2000	135	225	29/32	19/26	3
1632415041956	5000	180	316	45/40	24/29	1
1632415041957	5000	180	316	45/40	29/32*	1
1632415041966	10000	225	380	50/42	29/32	1
1632415041967	15000	260	420	50/42	29/32	1
1632415041968	20000	278	502	50/42	29/32	1

\* Outside the regulation

\* Abmessungen außerhalb der Norm

PN  
70 4606



## BOTTLES WOULFF, WITH GROUND

in SJ stopper and outlet stopcock

Bottles suitable for abstersion.

## ABKLÄRFLASCHEN

mit NS – Stopfen und NS – Hahn am Boden

Flaschen für Reinigung geeignet.

Cat. No. 2038				
Catalog No. Katalog-Nr.	Capacity (ml) Kapazität	h (mm)	SJ	Packing (pcs) Packung (Stück)
1632415046500	500	19/26	19/26	4
1632415046940	1000	24/29	19/26	4
1632415046950	2000	29/32	19/26	3
1632415046956	5000	45/40	24/29	3
1632415046966	10000	50/42	29/32	1
1632415046967	15000	50/42	29/32	1
1632415046968	20000	50/42	29/32	1

