

Información general

Nombre de producto: PICNÓMETRO GAY-LUSSAC CON TAPÓN

Fecha Rev.: Dic 12, 2023

Detalles:

Frasco de densidad calibrada, fabricado con vidrio de borosilicato, tipo Gay-Lussac

Calibrado para contener (TC, In)

Incluye certificado individual

Tapón NS 10/19 con capilar

Parte superior del tapón nivelada y pulida

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Calibración
0100521400	5	Sin calibrar

Ficha de producto

Ref.	Capacidad (ML)	Calibración
0100521401	10	Sin calibrar
0100521402	25	Sin calibrar
0100521403	50	Sin calibrar
0100521404	100	Sin calibrar
0100521500	5	Calibrado con certificado individual
0100521501	10	Calibrado con certificado individual
0100521502	25	Calibrado con certificado individual
0100521503	50	Calibrado con certificado individual
0100521504	100	Calibrado con certificado individual



General Information

Product name: PYCNOMETERS GAY-LUSSAC, WITH STOPPER

Rev. Date: Dec 12, 2023

Details:

Pyknometer, made of borosilicate glass, Gay-Lussac type

Calibrated to contain (TC, In)

Includes individual certificate

NS 10/19 stopper with capillary tube

Leveled and polished top of the stopper

Picture



Technical Information

Product datasheet

Ref.	Capacity (ML)	Calibration
0100521400	5	Unadjusted
0100521401	10	Unadjusted
0100521402	25	Unadjusted
0100521403	50	Unadjusted
0100521404	100	Unadjusted
0100521500	5	Adjusted with single certificate
0100521501	10	Adjusted with single certificate
0100521502	25	Adjusted with single certificate
0100521503	50	Adjusted with single certificate
0100521504	100	Adjusted with single certificate

