

Catalog of medical and laboratory products





MEDILAB ČR s.r.o. has been a leading Czech manufacturer and distributor of medical and laboratory consumables for over two decades.

With automated production processes, a strong commitment to product quality, and exceptional customer service, we are not only enhancing our market share in the Czech Republic but also expanding our reach internationally.



Index

Syringes

3-part syringes

Hypodermal needles

Infusion sets

Infusion set

Infusion set with port

Extension tubes

Extension tubes

2-way extension tubes

Transfusion set

Safety connectors and stopcocks

3-way stopcocks

Flow-in — safety injection site adapter

Urinal bags

Urinal bags

Urine bag holders

Urinal bag — leg type

Fastening stripes for urinal bags — leg type

Foley catheters

2-way Foley catheters Nelaton for men and women

2-way Foley catheters Tiemann

2-way Foley catheters Nelaton for women

2-way Foley catheters Nelaton for children

Closure for Foley catheter

Oxygen masks

Oxygen nasal cannulas and airways

Bedding and sheets

Disposable blanket

Disposable pillow

Disposable nonwoven sheets

Bed sheets

Gloves and medical brush

Latex examination gloves

Medical cleaning brush

Surgical face masks

Surgical caps

Surgical operating clothes

Surgical operating shirts

Surgical operating pants

Sampling brushes

Cytology brushes

Cervical brushes

Laboratory equipment

Microscope slides

Pipettes

Microtubes

Tubes

Centrifuge tubes

Microscope slide boxes

Cardboard slide mailers

Archiving box for histological cassettes

Container tray

Specimen containers

Specimen containers

Round bucket white

HistoMED

Pre-filled tubes HistoMED

Pre-filled containers HistoMED

Formalin

HistoMED

Surgical suture material

Syringes






Insulin syringes


The disposable insulin syringes come with an integrated 29G needle and are suitable for the administration of insulin.

3-part syringes

Disposable syringes are equipped with a double gasket, which prevents air penetration and ensures smooth movement of the plunger in the syringe cylinder. They are packed in a practical sterile peel-bag. The syringes are used for intravenous and intramuscular administration of medications by a needle.

Product description	REF	pcs/box	pcs/carton
 Insulin syringe U100 – 0,5 ml	IS0529G	100	3 600
 Insulin syringe U100 – 1 ml	IS1029G	100	3 600
 1 ml, sterile, centric cone, Luer-slip, non-residual	KD1301B	100	3 600
 2 ml, sterile, centric cone, Luer-slip	KD1302	100	2 700

3-part syringes





3-part syringes

Disposable syringes are equipped with a double gasket, which prevents air penetration and ensures smooth movement of the plunger in the syringe cylinder. They are packed in a practical sterile peel-bag.

Irrigation syringes

Disposable syringes are equipped with a double gasket, which prevents air penetration and ensures smooth movement of the plunger in the syringe cylinder. They are packed in a practical sterile peel-bag.



Product description	REF	pcs/box	pcs/carton
 50 ml, sterile, centric cone, Luer-lock	KD1352	50	300
 60 ml, sterile, centric cone, with cone	MA1450P	25	200

Hypodermal needles



Hypodermal needles

Disposable hypodermal needles are color-coded according to international codes. Needles are packed in sterile peel-bag.

Product description			REF	pcs/box	pcs/carton
 25G x 1"	0,50 mm x 25 mm	orange	HN25G	100	10 000
 23G x 1"	0,60 mm x 25 mm	blue	HN23G	100	10 000



Infusion sets




Infusion sets (length 150 cm)

Infusion sets are suitable for administering infusion solutions from bottles or plastic pouches into the patient's vein by gravity.

Infusion sets are packed in a single sterile package and are intended for single use. Without phthalates.




Product description	REF	pcs/box	pcs/carton
 Infusion set for gravity infusion	MA2208P	25	400
 Infusion set for gravity infusion with port	MA1108PORT	25	400

Extension tubes




Extension tubes

Extension tubes can be used in combination with infusion sets, transfusion sets and peripheral intravenous catheters. Without phthalates.


Product description	REF	pcs/box	pcs/carton
 Extension tubes 1,8 × 450	SH45	100	800
 Extension tubes 1,8 × 1800	SH180	50	400
 2-way extension tubes	SHY	100	800

Transfusion set



Transfusion set

Transfusion sets are particularly suitable for the intravenous transfusion of blood into the patient's vein. Without phthalates.

Product description	REF	pcs/box	pcs/carton
 Transfusion set	KD301N	25	250

Safety connectors and stopcocks





3-way stopcocks

Multi-way ventile for infusion sets and monitoring, lipidoresistant.

Flow-in – safety injection site adapter

Automatic valves are used for closing peripheral intravenous catheters and connection tubes. They facilitate drug administration.



Product description	REF	pcs/box	pcs/carton
 3-way stopcocks	STF11	50	1 000
 Flow-in	FLW	100	2 000

Urinal bags



Urinal bags

Urinal bags are sterile, with an antireflux valve, a capacity 2000 ml and a tube length of 90 or of 120 cm. The bag without outlet must be emptied after reaching its capacity and discarded in the determined place. The bag with outlet can be emptied during use using the cross valve. After use, it must be emptied and discarded in the determined place. All bag types have suspension holes for hanging the bag with the bag holder on the patient's bed.

Product description	REF	pcs/box	pcs/carton
 Urinal bags with bottom cross valve, 2000 ml, length of tube 90 cm, sterile single packing	UBT2000	10	250
 Urinal bags with bottom cross valve, 2000 ml, length of tube 120 cm, sterile single packing	UBT2000L	10	250
 Urinal bags without outlet, 2000 ml, length of tube 90 cm, single packing	UBT2000B	10	250
 Urine bag holder - hanger	UPHR	50	500
 Urine bag holder	UBH	100	2 000

Urinal bags





Urinal bag – leg type

Urinal bags have three chambers, an antireflux valve, capacity 750 ml and tube length of 45 cm. They comply with the needs of everyday usage. The bags are fastened with a reusable set of two velcro straps, which is not included in the packing. The rear part of the bag is textile-coated for sweat absorption and skin comfort. Urine is distributed in three chambers, keeping the bag flat for maximum discretion. The distribution of content in the chambers helps ensure secure fastening during movement.

Fastening stripes for urinal bags – leg type

Set of two velcro stripes for fastening the urinal bag to patient's leg. Stripes are coated with silicone stripe, which prevents slipping off. Fastening stripes are made of flexible latex-free material.

Product description	REF	pcs/box	pcs/carton
 Urinal bag – leg type, capacity 750 ml, three chambers, length of tube 45 cm, single packing	ULB500	10	200
 Set of fastening velcro stripes for urinal bags – leg type, length 50 cm, width 2,5 cm, white	LBS	10 pairs	500 pairs

Foley catheters




2-way Foley catheters Nelaton (long, length 40 cm, for men and women)

The permanent balloon catheter is sterile and 100% silicone coated. The integrated plastic valve allows easy inflation of the balloon.

Product description	REF	pcs/box	pcs/carton
○ Size CH12, balloon capacity 5 – 15 ml, for men and women	2FC12FRM5	10	500
● Size CH14, balloon capacity 5 – 15 ml, for men and women	2FC14FRM5	10	500
● Size CH16, balloon capacity 5 – 15 ml, for men and women	2FC16FRM5	10	500
● Size CH18, balloon capacity 5 – 15 ml, for men and women	2FC18FRM5	10	500
● Size CH20, balloon capacity 5 – 15 ml, for men and women	2FC20FRM5	10	500
● Size CH22, balloon capacity 5 – 15 ml, for men and women	2FC22FRM5	10	500
● Size CH24, balloon capacity 15 – 30 ml, for men and women	2FC24FRM5	10	500
● Size CH26, balloon capacity 15 – 30 ml, for men and women	2FC26FRM5	10	500

Foley catheters




2-way Foley catheters Tiemann (long, length 40 cm, only for men)

The permanent balloon catheter is sterile and 100% silicone coated. The integrated plastic valve allows easy inflation of the balloon.

Product description	REF	pcs/box	pcs/carton
○ Size CH12, balloon capacity 5 – 15 ml	2FC12TM	10	500
● Size CH14, balloon capacity 5 – 15 ml	2FC14TM	10	500
● Size CH16, balloon capacity 5 – 15 ml	2FC16TM	10	500
● Size CH18, balloon capacity 5 – 15 ml	2FC18TM	10	500
● Size CH20, balloon capacity 5 – 15 ml	2FC20TM	10	500
● Size CH22, balloon capacity 5 – 15 ml	2FC22TM	10	500
● Size CH24, balloon capacity 5 – 15 ml	2FC24TM	10	500

Foley catheters




2-way Foley catheters Nelaton for women (length 20 cm, only for women)

The permanent balloon catheter is sterile and 100% silicone coated. The integrated plastic valve allows easy inflation of the balloon.

Product description	REF	pcs/box	pcs/carton
○ Size CH12, balloon capacity 5 – 10 ml, for women	2FC12FRW	10	500
● Size CH14, balloon capacity 5 – 10 ml, for women	2FC14FRW	10	500
● Size CH16, balloon capacity 5 – 15 ml, for women	2FC16FRW	10	500
● Size CH18, balloon capacity 5 – 15 ml, for women	2FC18FRW	10	500
● Size CH20, balloon capacity 5 – 15 ml, for women	2FC20FRW	10	500
● Size CH22, balloon capacity 5 – 15 ml, for women	2FC22FRW	10	500

Foley catheters



2-way Foley catheters Nelaton for children (length 20 cm)

The permanent balloon catheter is sterile and 100% silicone coated. The integrated plastic valve allows easy inflation of the balloon.

Closure for Foley catheter „spigot“

Intended for closing of foley catheter (Nelaton, Tiemann).






Product description	REF	pcs/box	pcs/carton
● Size CH8, baloon capacity 3 ml, for children	2FC8FRM3	10	500
● Size CH10, baloon capacity 3 ml, for children	2FC10FRM3	10	500
● Closure for Foley catheter - „spigot“, sterile	UFK	100	1 000

Oxygen masks



Oxygen masks — standard / nonrebreathable / with nebulizer

Oxygen masks are sterile, with a tube length of 213 cm.

Product description	REF	pcs/package	pcs/carton
 Oxygen mask — standard — size L	KMSL	1	100
 Oxygen mask — standard — size XL	KMSXL	1	100
 Oxygen masks — nonrebreathable, size L	KMVL	1	50
 Oxygen masks — nonrebreathable, size XL	KMVXL	1	50
 Oxygen mask — nebulizer — size L	KMRL	1	100

Oxygen nasal cannulas and airways



Oxygen nasal cannulas

Oxygen nasal cannulas for adults. The tube is equipped with cross-section, which prevents it from its complete cranking. Length 213 cm. Without phthalates.

Airways

Airways are anatomically shaped, compact, disposable and sterile. Different colors for size distinction.

Product description	REF	pcs/package	pcs/carton
 Oxygen nasal cannulas – for adults	KB	20	200
 Airways – size 3, length 90 mm, yellow	GA3	1	50
 Airways – size 4, length 100 mm, red	GA4	1	50

Bedding and sheets



Disposable blanket

Disposable blanket with polyester filling is made of non-woven fabric. The blanket is waterproof.

Disposable pillow

The disposable pillow with polyester filling is made of non-woven fabric. It is particularly suitable for ARO, ICU, infection departments or rescue services.

Disposable nonwoven sheets

The sheets are made of non-woven fabric. They are intended to use as a one-time hygienic protection for hospital bedding.

Disposable transport sheet

The disposable transport sheet with 12 stamped-in carrying handles is made of strong polyester fabric. It is used for transporting immobilized patients and transferring them from one surface to another. Suitable for patients up to 150 kg.

Product description	REF	pcs/box	pcs/carton
● Disposable blanket 190 × 110 cm, 300 g, blue	DJ300B	1	54
● Disposable blanket 190 × 110 cm, 400 g, blue	DJ400B	1	48
● Disposable blanket 190 × 110 cm, 500 g, blue	DJ500B	1	40
● Disposable nonwoven blanket cover 195 × 115, 25 g	JP25195Z	10	100
● Disposable pillow 40 × 40 cm, white	PJ40x40-W	1	36
● Disposable pillow-case 42 × 42 cm, 25 g	JP2542	10	200
● Disposable transport sheet 200 × 70 cm, 230 g/m ²	RTETT	1	30

Bed sheets



Disposable bed sheets - nonwoven textile

Bed sheets can be used as disposable hygienic protection on hospital and examination beds. Bedsheet REF PRO3015200 is additionally equipped with polyethylene waterproof layer.

Product description	REF	pcs/package	pcs/carton
Disposable bed sheet, nonwoven textile, 220 × 80 cm, 25 g/m ²	PRO25220	10	100
Disposable bed sheet, nonwoven textile, 220 × 120 cm, 25 g/m ²	PRO25220120	10	100
Disposable bed sheet, nonwoven textile, 220 × 140 cm, 25 g/m ²	PRO25240140	10	100
Disposable bed sheet, nonwoven textile, 200 × 80 cm, PP 30 g/m ² + PE 15 g/m ²	PRO3015200	10	100
Disposable bed sheet, nonwoven textile, 200 × 80 cm, 40 g/m ²	PRO40200	10	100
Disposable bed sheet, nonwoven textile, 220 × 140 cm, 40 g/m ²	PRO40220140	10	100
Disposable bed sheet, nonwoven textile, 200 × 90 cm, 60 g/m ²	PRO60200	10	100
Disposable bed sheet, nonwoven textile, 210 × 90 cm, 60 g/m ²	PRO60210	10	100
Disposable bed sheet, nonwoven textile, 220 × 90 cm, 60 g/m ²	PRO60220	10	100
Disposable bed sheet, nonwoven textile, 220 × 140 cm, 60 g/m ²	PRO60220140	5	50

Gloves and medical brush



Latex examination gloves

Examination gloves are left-right and they are for single use. They distinguish in high durability as well as in required softness. They are used for examination, nursing, diagnostic, terapeutical and other procedures in health sector.

Medical brushes

The medical brush has soft bristles, which reliably clean joints, nails and skin while being very gentle. The brush has integrated sponge, which reliably cleans hands and forearms without causing irritation or damage (microtrauma) to the skin. The package also includes a nail cleaner.

Product description	REF	pcs/box	pcs/carton
 Latex examination gloves — size S, powdered	LEGSPP	100	1 000
 Latex examination gloves — size M, powdered	LEGMPP	100	1 000
 Latex examination gloves — size L, powdered	LEGLPP	100	1 000
 Latex examination gloves — size S, powder-free	LEGSPPF	100	1 000
 Latex examination gloves — size M, powder-free	LEGMPF	100	1 000
 Latex examination gloves — size L, powder-free	LEGLPF	100	1 000
 Medical cleaning brush with sponge and nail cleaner, sterile	MB	30	300

Surgical face masks



Surgical face masks



The bacterial efficiency of surgical face masks is $\geq 98\%$.

The material is non-irritating and comfortable for the skin.

Face masks with ties are fastened with binding ties and face mask with elastic strips is fastened with comfortable elastic strips.

Both types of face masks are equipped with an adjustable strip in the nose and mouth area.

All face masks comply with standard EN 14683. Protection type II..

Product description	REF	pcs/box	pcs/carton
 3-ply face mask with ties, blue	CHOM1	50	2 000
 3-ply face mask with elastic strips, blue	CHOM2	50	2 000

Surgical face masks



Surgical face masks — with flowers / with hearts / black with silver flowers

The bacterial efficiency of surgical face masks is $\geq 98\%$.

The material is non-irritating and comfortable for the skin.

Face masks with ties are fastened with binding ties and face mask with elastic strips is fastened with comfortable elastic strips.

Both types of face masks are equipped with an adjustable strip in the nose and mouth area.

All face masks comply with standard EN 14683. Protection type II.

Product description	REF	pcs/box	pcs/carton
 3-ply face mask with ties, blue flowers	CHOM5	50	2 000
 3-ply face mask with ear-loops, blue flowers	CHOM6	50	2 000
 3-ply face mask with ties, pink flowers	CHOM7	50	2 000
 3-ply face mask with ear-loops, pink flowers	CHOM8	50	2 000
 3-ply face mask with ties, pink hearts	CHOM9	50	2 000
 3-ply face mask with ear-loops, pink hearts	CHOM10	50	2 000
 3-ply face mask with ties, black with silver flowers	CHOM11	50	2 000
 3-ply face mask with ear-loops, black with silver flowers	CHOM12	50	2 000

Surgical caps



Surgical cap – with ties




Disposable blue surgical cap with ties is made of soft breathable nonwoven textile. Surgical cap is suitable for using during surgical operation as well as universal head cover.

Surgical cap – with ties and absorption strip

Disposable blue surgical cap with ties is made of soft breathable nonwoven textile. Surgical cap is equipped with absorption stripe. It is suitable for using during surgical operation as well as universal head cover.

Surgical cap – with elastic strip

Disposable blue surgical cap is made of soft breathable nonwoven textile. Surgical cap has elastic strip in the rear part and is suitable for using during surgical operation as well as universal head cover.

Product description	REF	pcs/box	pcs/carton
 Surgical cap with ties, blue	SCTOB	100	1 000
 Surgical cap with ties and absorption strip, blue	SCWAS	100	1 000
 Surgical cap with elastic strip, blue	SCPB	100	1 000

Surgical caps



Surgical caps – with flowers / with hearts / black with silver flowers

Disposable surgical caps are available with ties or elastic strip. They are made of soft breathable nonwoven textile. Surgical caps are suitable for using during surgical operation as well as universal head cover.

Product description	REF	pcs/box	pcs/carton
 Surgical cap with ties, blue flowers	SCTOPWB	100	1 000
 Surgical cap with ties, pink flowers	SCTOPWP	100	1 000
 Surgical cap with ties, pink hearts	SCTOPH	100	1 000
 Surgical cap with ties, black with silver flowers	SCTOPBS	100	1 000
 Surgical cap with elastic strip, blue flowers	SCPWB	100	1 000
 Surgical cap with elastic strip, pink flowers	SCPWP	100	1 000
 Surgical cap with elastic strip, pink hearts	SCPH	100	1 000
 Surgical cap with elastic strip, black with silver flowers	SCPBS	100	1 000

Surgical operating clothes



















Surgical operating shirt — flowered pocket / pocket with hearts


Disposable surgical shirts are made of soft, breathable non-woven fabric. Shirts are suitable to be worn during surgery. They are available in sizes S, M, L, XL and variants with a flowered pocket or a pocket with hearts. The shirts have short sleeves and a V-neck.

Surgical operating pants

Disposable surgical pants are straight-cut with a tie at the waist, are made of soft, breathable non-woven fabric. Pants are suitable to be worn during surgery and available in sizes S, M, L and XL.

Product description	REF	pcs/box	pcs/carton
 Surgical operating shirt — pink flowers S	OSSPF		50
 Surgical operating shirt — pink flowers M	OSMPF		50
 Surgical operating shirt — pink flowers L	OSLPF		50
 Surgical operating shirt — pink flowers XL	OSXLPF		50
 Surgical operating shirt — blue flowers S	OSSBF		50
 Surgical operating shirt — blue flowers M	OSLBF		50
 Surgical operating shirt — blue flowers L	OSMBF		50
 Surgical operating shirt — blue flowers XL	OSXLBF		50
 Surgical operating shirt — pink hearts S	OSSH		50
 Surgical operating shirt — pink hearts M	OSMH		50
 Surgical operating shirt — pink hearts L	OSLH		50
 Surgical operating shirt — pink hearts XL	OSXLH		50
 Surgical operating pants — blue S	OTS		50
 Surgical operating pants — blue M	OTM		50
 Surgical operating pants — blue L	OTL		50
 Surgical operating pants — blue XL	OTXL		50

Sampling brushes







Cytology brushes

Soft sampling brush for taking cytological samples.

Cervical brushes

Soft sampling brush used to collect samples from the cervix.

Product description	REF	pcs/box	pcs/carton
 Cytology brush — sterile	CBS	100	5 000
 Cytology brush — nonsterile	CB	100	5 000
 Cytology brush — nonsterile, with ball tip	CBK	100	5 000
 Cytology brush — nonsterile, with big head and ball tip	CBKT	100	5 000
 Cytology brush — cervical	CBCR	100	2 500
 Cervical brush — hardness 12 %, nonsterile	CBCH2	100	2 500
 Cervical brush — hardness 15 %, nonsterile	CBCH3	100	2 500

Laboratory equipment



Microscope cover glass

The slides are made of high-quality clear glass and are vacuum packed by 500 pieces.

Microscope slides — rounded on the both sides

Clear or one side frosted, resistant to commonly used laboratory solvents.

Product description	REF	pcs/box	pcs/carton
 Microscope cover glass — corners 90°, dimensions 24 × 50 mm, thickness 0,13 — 0,17 mm	MCG7201	100	20 000
 Microscope slides, rounded on the both sides, pre-cleaned — clear, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7101	50	2 500
 Microscope slides, rounded on the both sides, one side frosted, pre-cleaned — corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7105-1	50	2 500

Laboratory equipment









Microscope slides — printable


These slides are printable and rounded on the both sides.

Microscope slides — printable, positive charged

These slides have a permanent positive charge. The static electricity attracts parts of frozen tissues and cytology preparations onto the slide surface. Parts of tissues and cytology preparations adhere better to positive charged slides, without the need of special glue or protein coating. Positively charged slides also include adhesive microscope slides that are coated with silane.

Product description	REF	pcs/box	pcs/carton
 Microscope slides, printable, pre-cleaned — white, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7109-W	50	2 500
 Microscope slides, printable, pre-cleaned — white, corners 90°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	ME7109-W90	50	2 500
 Microscope slides, positive charged, printable, pre-cleaned — green, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEPN-G	50	2 500
 Microscope slides, positive charged, printable, pre-cleaned — pink, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEPN-P	50	2 500
 Microscope slides, positive charged, printable, pre-cleaned — yellow, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEPN-Y	50	2 500
 Adhesive microscope slides, printable, coated with silane, positive charged, pre-cleaned — blue, corners 45°, dimensions 25,4 × 76,2 mm, thickness 1,0 — 1,2 mm	MEAS-B	50	2 500

Laboratory equipment



Pipettes

Disposable pipette is made of special nontoxic low density polyethylene. It is suitable for safe transfer and distribution of liquids in all laboratory types. Eliminated risk of infection transfer.

Microtubes









Microtubes are made of robust virgin polypropylene and are transparent for better visibility. Tubes with a safety lock have perfect stability, tightness, protection against vaporization and unwanted opening. The flat lid with tapered edges enables easy opening and closing. The frosted part on the side is intended for easy description or marking, and the side is also equipped with a scale. Microtubes are autoclavable and microtubes code 35200N are suitable for use in a centrifuge.

Tubes

The tubes provide perfect tightness and protection against evaporation and accidental opening. Red or blue polyethylene stopper is available for the test tubes.

Centrifuge tubes

The test tubes are intended for the preparation, processing and storage of biological samples. The tubes are equipped with a scale and a surface for easy labeling. The screw cap ensures a perfect seal. The tubes are autoclavable and intended for centrifugation.

Product description	REF	pcs/package	pcs/carton
 Disposable pipette, 4 ml	13301	500	3 000
 Disposable graduated pipette, 1 ml	JPG01	500	5 000
 Microtube – 1,5 ml, without lid	36210	1 000	11 000
 Microtube with scale, with lid	35200N	500	10 000
 Test tube 16 x 80 mm, capacity 10 ml, freestanding, with lid	86120	100	1 000
 Tube 16 x 100 mm, volume 14 ml, without lid	81621	500	2 500
 Polyethylene lid for tube, red color	89216C	1 000	16 000
 Polyethylene lid for tube, blue color	89216M	1 000	16 000
 Tube with conical bottom, free standing, volume 50ml, with lid	HP1001	25	500
 Tube with conical bottom, volume 50ml, with lid	HP1003	25	500

Laboratory equipment



Microscope slide boxes

Plastic boxes are used to store and transfer microscope slides and prevent them from breaking or sticking.

Cardboard slide mailers











Slide mailers are made of cardboard and designed for holding microscope slides.

Archiving box for histological cassettes

Cardboard box is designed for storing histological cassettes.

Container tray

The polystyrene tray is used to store containers with a diameter of 45 and 32 mm.

Product description	REF	pcs/package	pcs/carton
 Microscope slide box for 5 pcs of slides, transparent	HP4007	4	100
 Microscope slide box for 12 pcs of slides (orange colour)	HP4008	2	
 Microscope slide box for 50 pcs of slides	BFCS	1	100
 Microscope slide box for 100 pcs of slides	BFCS-L	1	50
 Slide mailer 295 x 205 mm	DEKLS	10	40
 Slide mailer 230 x 460 mm	DEKLM	10	30
 Slide mailer 230 x 485 mm	DEKLL	10	30
 Slide mailer — pad without lid	DEKL-BV	10	40
 Archiving box for histological cassettes 300 x 420 x 50 mm	AKE300x420x50		50
 Polystyrene tray for containers with diameter 45 and 32 mm	XPS255x205x50	1	

Specimen containers



Specimen containers

Specimen containers are made of nontoxic polypropylene. The lids are made of HDPE or polypropylene. The containers are suitable for collection, transport and storage of samples and prevent leakage and evaporation of liquids. They prevent leakage and vaporization of formalin. Upper and bottom grooves ensure easy opening and closing with gloves. The scales are placed on three sides of containers.

Product description	REF	pcs/carton
 Specimen containers, capacity 20 ml	28100	1 000
 Specimen containers, capacity 40 ml	28200	500
 Specimen containers, capacity 60 ml	28300	500
 Specimen containers, capacity 90 ml	28400	300
 Specimen containers, capacity 120 ml	28500	300
 Specimen containers, capacity 160 ml	28900	200
 Specimen containers, capacity 250 ml	28600	154
 Specimen containers, capacity 500 ml	28700	132
 Specimen containers, capacity 1000 ml	28800	64
 Specimen containers, capacity 7 ml	22500	700
 Specimen containers, capacity 30 ml	22100	400

Specimen containers



Round bucket white

Buckets with lids are made of polypropylene and are suitable for collecting, transporting and storing samples.

Product description	REF	pcs
 Round bucket white 1 l, plastic handle + lid	EZ103	1
 Round bucket white 3 l, plastic handle + lid	EZ301	1
 Round bucket white 5 l, plastic handle + lid	EZ504	1
 Round bucket white 8,5 l, plastic handle + lid	EZ850	1
 Round bucket white 10 l, plastic handle + lid	EZ1002	1
 Round bucket white 20 l, metal handle + lid with safety lock	EZ2002	1

HistoMED



Pre-filled tubes and containers HistoMED

HistoMED is a sample container always half-filled with a fixing solution intended for the fixation, storage and transport of biological material. The container contains 10% buffered stabilized formalin with a formaldehyde content of max. 4%, pH 7.2 – 7.4.

*) minimum order quantity

Product description	REF	pcs/package	pcs/carton
 Pre-filled tube HistoMED, volume 1,5 ml, amount of formalin 0,7 ml	HMED05		1 000*
 Pre-filled tube HistoMED, volume 7 ml, amount of formalin 3,5 ml	HMED07	50*	1 000
 Pre-filled tube HistoMED, volume 14 ml, amount of formalin 7 ml	HMED14	100	800*
 Pre-filled container HistoMED, vol. 20 ml, amount of formalin 10 ml	HMED20	60	120*
 Pre-filled container HistoMED, vol. 40 ml, amount of formalin 20 ml	HMED40	35	140*
 Pre-filled container HistoMED, vol. 60 ml, amount of formalin 30 ml	HMED60	35	105*
 Pre-filled container HistoMED, vol. 90 ml, amount of formalin 45 ml	HMED90	35	70*
 Pre-filled container HistoMED, vol. 120 ml, amount of formalin 60 ml	HMED120	23	46*
 Pre-filled container HistoMED, vol. 160 ml, amount of formalin 80 ml	HMED160	23	23*
 Pre-filled container HistoMED, vol. 250 ml, amt. of formalin 125 ml	HMED250	9	27*
 Pre-filled container HistoMED, vol. 500 ml, amt. of formalin 250 ml	HMED500		10*
 Pre-filled container HistoMED, vol. 1000 ml, amt. of formalin 500 ml	HMED1000		5*
 Pre-filled container HistoMED, vol. 3000 ml, amt. of formalin 1500 ml	HMED3000	1*	8
 Pre-filled container HistoMED, vol. 5000 ml, amt. of formalin 2500 ml	HMED5000	1*	4

HistoMED








Formalin

The bottle or canister is completely filled with 10% buffered stabilized formalin with a maximum formaldehyde content of 4%, pH 7.2 - 7.4.

HistoMED

HistoMED is a sample container always half filled with fixing solution for fixation, storage and transport of biological material. The container contains 32% stabilized formalin with a formaldehyde content of max. 10%, pH 9.5 - 10.5.

Product description	REF	pcs/package	pcs/carton
 Buffered formalin in a 1000 ml bottle	FORMAL1	1	8
 Buffered formalin in 5000 ml canister	FORMAL5	1	3
 HistoMED prefilled container, volume 500 ml, amount of formalin 250 ml	HMED10500		10
 HistoMED prefilled container, volume 1000 ml, amount of formalin 500 ml	HMED101000		5
 HistoMED prefilled container, volume 3000 ml, amount of formalin 1500 ml	HMED103000	1	8
 HistoMED prefilled container, volume 5000 ml, amount of formalin 2500 ml	HMED105000	1	4

Absorbable surgical suture material



	<p>Surgicryl® 910 Polyglactine 910</p> <p>Support of wound approximately 30 days</p>		
	<p>Surgicryl® PGA Polyglycolic acid</p> <p>Support of wound approximately 30 days</p>		
	<p>Surgicryl® Rapid Polyglycolic acid</p> <p>Support of wound approximately 10 days</p>		
	<p>Surgicryl® Monofast Polyglecaprone 25</p> <p>Support of wound approximately 20 days</p>		
	<p>Surgicryl® Monofilament Polydioxanone</p> <p>Support of wound approximately 60 days</p>		

Non-absorbable surgical suture material



You can request the catalog from our sales representatives or order by e-mail objednavky@medilab.cz. For more information, visit www.medilab.cz.



Quality Equipment for Medical

Certification



Certifikát ES
Směrnice 93/42/EHS Příloha V
Zajištění jakosti výroby
zdravotnických prostředků

Registrační číslo: DD 60144023 0001
Zpráva z auditu číslo: 26300449 002

Výrobce: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká republika

Výrobky: (zařazené výrobky jsou uvedeny v příloze)

Datum vypršení platnosti: 2024-04-08


Notifikovaný orgán:
D. Swiatko

Datum vydání: 2019-11-22

Datum vystavení: 2019-11-22

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg

TÜV Rheinland LGA Products GmbH je Notifikovaný orgán podle směrnice 93/42/EHS o zdravotnických prostředcích. Jako identifikátor číslo Notifikovaného orgánu je 0197.



TÜV Rheinland
LGA Products GmbH
Tillystraße 2, 90431 Nürnberg

Doc. 1/3, Rev. 0

Příloha k certifikátu

s registračním číslem: DD 60144023 0001
Zpráva z auditu č.: 26300449 002

Výrobce: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká republika

Zařazené výrobky:

- Sterilní nasazovací ústevky Flow-in
- Sterilní 3-cestné uzavírací kobotky
- Sterilní Foleyovy katetry HELATOR
- Sterilní infuzní stříkačky
- Sterilní infuzní sety
- Sterilní nasazovací ústevky
- Sterilní nasazovací ústevky
- Sterilní Foleyovy katetry SIMON

Pro následující zdravotnické prostředky pokrývá certifikace pouze aspekty výroby spojené se zajištěním a udržováním sterility prostředků:

- Chirurgické kartáčky
- Třísnivé sáčky

Notifikovaný orgán:
D. Swiatko

Datum: 2019-11-22



Certificate
Quality Management System
EN ISO 13485:2016

Registration No.: SX 1550678-1
Organization: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Czech Republic

Scope: Production and distribution of sterile and non-sterile disposable medical devices.

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the abovementioned standard are fulfilled. The quality management system is subject to periodic surveillance.

Report No.: 8490384-20
Effective date: 2022-09-13
Expiry date: 2025-04-08
Issue date: 2022-09-13

TÜV Rheinland LGA Products GmbH
Tillystraße 2 - 90431 Nürnberg, Germany

1/1



EC Certificate
Directive 93/42/EEC Annex V
Production Quality Assurance
Medical Devices

Registration No.: DD 60144023 0001
Report No.: 26300449 002

Manufacturer: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká republika


Products: (see attachment for products included)

Expiry Date: 2024-04-08

Notified Body:
D. Swiatko

Datum: 2019-11-22

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



TÜV Rheinland
LGA Products GmbH
Tillystraße 2, 90431 Nürnberg

Doc. 1/3, Rev. 0

Attachment to Certificate
Registration No.: DD 60144023 0001
Report No.: 26300449 002

Manufacturer: **MEDILAB ČR s.r.o.**
Záhumná 574
51771 Česká Meziříčí
Česká republika

Products included:

- Sterile injection site adapters Flow-in
- Sterile 3-way stopcocks
- Sterile Foley catheters HELATOR
- Sterile insulin syringes
- Sterile syringes for single use
- Sterile infusion sets
- Sterile transfusion sets
- Sterile extension tubes
- Sterile Foley catheter SIMON

For the following medical devices the scope covers only the aspects of manufacture concerned with securing and maintaining sterile conditions:

- Medical brushes
- Drain bags

Notified Body:
D. Swiatko

Date: 2019-11-22

Certificates of quality management system ČSN EN ISO 13485:2016.

Certificates of product quality according to 93/42 EEC, justification for product CE marking.



TÜVRheinland
CERTIFIED

Management System
EN ISO 13485:2016
+AC:2018+A11:2021



www.tuv.com
ID 9000023503



MEDILAB®

MEDILAB ČR s.r.o., Czech Republic
E-mail: medilab@medilab.cz, objednavky@medilab.cz
Phone: (+420) 494 661 067, Fax: (+420) 494 661 068

www.medilab.cz