


Pre-filled containers for biological samples

HistoMED



Quality Equipment for Medical


MEDILAB ČR s.r.o., Záhumnská 574, 51771 České Meziříčí, Czech Republic
Phone: (+420) 494 661 067, Mobile: (+420) 774 711 573, E-mail: medilab@medilab.cz



Pre-filled containers / HistoMED

Pre-filled tubes and containers HistoMED

HistoMED is a sample container always half filled with fixing solution for fixation, storage and transport of biological material. The container contains 10% buffered stabilized formalin with a formaldehyde content of max. 4%, pH 7.2 - 7.4.

Product description	REF	pcs/package	pcs/carton
 Pre-filled tube HistoMED, volume 7 ml, amount of formalin 3,5 ml	HMED07	50	1 000
 Pre-filled tube HistoMED, volume 14 ml, amount of formalin 7 ml, coloured with eosin	HMED14	100	800
 Pre-filled container HistoMED, volume 20 ml, amount of formalin 10 ml	HMED20	60	120
 Pre-filled container HistoMED, volume 40 ml, amount of formalin 20 ml	HMED40	35	140
 Pre-filled container HistoMED, volume 60 ml, amount of formalin 30 ml	HMED60	35	105
 Pre-filled container HistoMED, volume 90 ml, amount of formalin 45 ml	HMED90	35	70
 Pre-filled container HistoMED, volume 120 ml, amount of formalin 60 ml	HMED120	23	46
 Pre-filled container HistoMED, volume 160 ml, amount of formalin 80 ml	HMED160	23	23
 Pre-filled container HistoMED, volume 250 ml, amount of formalin 125 ml	HMED250	9	27
 Pre-filled container HistoMED, volume 500 ml, amount of formalin 250 ml	HMED500		10
 Pre-filled container HistoMED, volume 1000 ml, amount of formalin 500 ml	HMED1000		5
 Pre-filled container HistoMED, volume 3000 ml, amount of formalin 1500 ml	HMED3000	1	8
 Pre-filled container HistoMED, volume 5000 ml, amount of formalin 3000 ml	HMED5000	1	4

Formalin

The bottle or canister is completely filled with 10% buffered stabilized formalin with a maximum formaldehyde content of 4%, pH 7.2 - 7.4.

HistoMED

HistoMED is a sample container always half filled with fixing solution for fixation, storage and transport of biological material. The container contains 32% stabilized formalin with a formaldehyde content of max. 10%, pH 9.5 - 10.5.

Product description	REF	pcs/package	pcs/carton
 Buffered formalin in a 1000 ml bottle	FORMAL1	1	8
 Buffered formalin in 5000 ml canister	FORMAL5	1	3
 Pre-filled container HistoMED, volume 500 ml, amount of formalin 250 ml	HMED10500	1	10
 Pre-filled container HistoMED, volume 1000 ml, amount of formalin 500 ml	HMED101000	1	5
 Pre-filled container HistoMED, volume 3000 ml, amount of formalin 1500 ml	HMED103000	1	8
 Pre-filled container HistoMED, volume 5000 ml, amount of formalin 2500 ml	HMED105000	1	4